YC-501 LABORATORY SCALE SPRAY DRYER

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
Since 2005, Shanghai Pilotech has provided mini spray dryer YC-015, laboratory low-temperature spray dryer YC-1800, pilot scale spray dryer YC-018, laboratory vacuum spray dryer YC-2000, laboratory inert loop spray dryer YC-015A and pilot scale inert loop spray dryer YC-018A system, having been the leader of laboratory scale spray dryer field in global market. We have installed over 3000 spray dryers in colleges, institutes and business R&D departments all over the world. In 2016, in accordance with the principle of providing reliable and durable products with accuracy and innovation, focusing on practical laboratory requirements of customers, Shanghai Pilotech launched brand-new YC-500 mini spray dryer, which was generated from many years of spray dryer manufacturing experience and modern instrument manufacturing. Thus, the user experience, the

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Spray drying can be used in a wide range of applications where the production of a free-flowing powder sample is required. This technique has successfully processed materials in the following areas:

Beverages • Flavours and Colourings • Milk and Egg Products • Plant and Vegetable Extracts • Pharmaceuticals • Heat Sensitive Materials • Plastics • Polymers and Resins • Perfumes • Ceramics and Advanced Materials • Soaps and Detergents • Blood • Dyestuffs • Foodstuffs • Adhesives • Oxides • Textiles • Bones, Teeth and Tooth Amalgam and many others. Most solutions and suspensions can be spray dried providing that the resulting product has the characteristics of a solid material.

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

A menu driven microprocessor controller allows the selection of inlet temperature, airflow, automatic de-blocker frequency and pump speed. The controller features an RS 232 output for connection to a PC or datalogger and software allows the control and monitoring of all functions and printing of results.

The self-priming peristaltic pump delivers the sample liquid from a container through a small diameter jet into the main chamber. At the same time an integral compressor pumps air into the outer tube of the jet which causes the liquid to emerge as a fine atomised spray into the drying chamber.

Heated air is blown through the main chamber evaporating the liquid content of the atomised spray. The solid particles of the material, which are normally in a free flowing state, are then separated from the exhaust air flow by a cyclone and collected in the sample collection bottle. The exhaust airflow is directed through a flexible 60 mm diameter hose direct to atmosphere or to an existing extraction system.

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TECHNICAL INFORMATION

全自动控制与手动控制双重控制模式，整个实验过程彩色触摸屏动态显示（动画）
Two modes of run: Automatic-mode, and Eye-monitored mode for the purpose of easily controlling experimental process.

进风温度控制：60℃ ~200℃
Air inlet temperature range 60°C to 200°C

出风温度控制：45℃ ~100℃
Air outlet temperature range 45°C to 100°C

物料处理量：最大500ml/h
Evaporation rate of water at inlet temperature of 200°C using
Standard: Approximately 500 ml/hour

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彩色LCD触摸屏参数显示：进风口温度/出风口温度/蠕动泵转速 /风量/通针频率
As standard a touch screen to display temperature, airflow volume, air pressure, pump speed and de-blocker frequency.

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Homocentric spray head: during experiments, the new
spray head would never be off-center in atomization,
and there won't be affected recovery rate by
spraying on the bottle wall. The spray head can
change the spray site to improve the spray effect by
moving up and down; notably, the spray head has
excellent adaptability to material with high
viscosity and few blocking situation;

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The first on-table model: the instrument can be placed on the experiment table, or Pilotech inert loop system YC-501, featuring small size and convenient operation;

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For aqueous solution and organic solvent: YC-500 mini spray dryer has the basic type and the organic solvent type YC-501 (with inert loop system). YC-500 is for aqueous solution, which can simply be updated to organic solvent type after equipped with YC-501 inert loop system;

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To upgrade the basic type into organic solvent type, YC-501 suit is the only thing to be brought: due to limited budget, many researchers with the requirement of organic solvent drying can only afford the aqueous solution type. However, they may later need to buy aqueous solution, which greatly wastes the scientific research funds. The YC-500 type laboratory on-table spray dryer completely solved this problem. The researchers only need to buy the YC-501 nitrogen cycle system to form the enclosed spray dryer after having the basic type. The connection is convenient. Using compressed air with aqueous solution and nitrogen for organic solvent, researchers will not be annoyed by these problems;

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Homocentric spray head: during experiments, the new spray head would never be off center in atomization, and there won't be affected recovery rate by spraying on the bottle wall. The spray head can change the spray site to improve the spray effect by moving up and down; notably, the spray head has excellent adaptability to material with high viscosity and few blocking situation;

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Tail gas dedusting equipment: YC-500 standard tail gas dedusting equipment can reduce environmental pollution while effectively improving the yield, especially the yield of difficult recycle;

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New type of heater: It is a challenge that laboratory spray dryer's heater strip can be easily burnout. The YC-500 effectively solved this problem and the researchers will not worry about changing the heater frequently;

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电子信箱：sales@pilotplant.cn

Business Center: Rm 1603, Block B, Everwin Building 521 South Wan Ping Rd., Shanghai, 200032, PR China

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SHANGHAI PILOTECH INSTRUMENT & EQUIPMENT CO., LTD